rev. 09/18

For research use only

## **Recombinant Protein A**

**CATALOG #:** 6500B-10 10 mg

6500B-25 25 mg 6500B-100 100 mg 6500B-1000 1 g

SOURCE: E. coli

**SEQUENCE:** Amino acid 32-327 of the Staphylococcus aureus subsp.

aureus. Protein A Ig binding domains with 6x His-tag on N-terminus. Gene Bank Accession Number YP\_498670.

**PURITY:** >97% by SDS-PAGE and HPLC analyses

MOL. WEIGHT: 35.5 kDa (Apparent M.W. 79 kDa (SEC))

THEORETICAL PI: 5.56

**EXTINCTION COEFFICIENT:** 7450 (at 280 nm)

FORM: Lyophilized

**RECONSTITUTION:** Reconstitute in water

**STORAGE CONDITIONS:** -20°C or below. Stable, as supplied, for at least 1 year.


hlgG BINDING: >95%

**BIOBURDEN:** No organisms detected

## APPLICATION:

The recombinant Protein A is a genetically engineered protein containing five Ig-binding regions of protein A and 6x His-tag on N-terminus. Cell wall binding region, albumin binding region and other non-specific binding regions have been eliminated from the recombinant Protein A to ensure the maximum specific IgG binding. The 6x His-tag can be used for affinity purification or for protein A detection using anti-His-tag antibody. The recombinant Protein A is ideal for purification of polyclonal or monoclonal IgG antibodies.

Protein A binds to most human and mouse IgG subclasses (e.g., human IgG1, IgG2, IgG4; mouse IgG1, IgG2a, IgG2b, IgG3). It also binds to total IgG from cow, guinea pig, hamster, horse, pig, and rabbit. Protein A has little affinity to chicken, goat, rat and sheep.



A: SDS-PAGE (4-20%) of Recombinant Protein A. Lane 1-3: 5, 2.5, and 1 μg, Protein A under Reducing conditions, Lane4: protein Ladder; Lane5-7: 5, 2.5, and 1 μg of Protein A under Non-Reducing conditions. **B: SEC chromatography of Protein A.** The apparent molecular weight is ~ 70 KDa Calculated Based on SEC standard curve (data not shown)

## **RELATED PRODUCTS:**

- Protein A (Lyophilized) (Cat. No. 6500B-10, -25, -100, -1000)
- Protein A-Agarose (Cat. No. 6526-1, -5, -25, -100)
- Hi-Bind™ Protein A-Agarose (Cat. No. 6520-1, -5, -25, -100)
- Protein A-Sepharose (Cat. No. 6501-1, -5, -25, -100)
- Hi-Bind™ Protein G-Agarose (Cat. No. 6513-1, -5, -25, -100)
- Protein G-Sepharose (Cat. No. 6511-1, -5, -25, -100)
- Protein A/G-Sepharose (Cat. No. 6503-1, -5, -25, -100)
- Protein A/G/L-Sepharose (Cat. No. 6541-1, -5, -25, -100)
- Protein A/G/L (Cat. No. 6540-1, -25, -100, -1000)

FOR RESEARCH USE ONLY! Not to be used on humans.