



Minute™ Lysosome Isolation Kit

(For mammalian cells/tissues)

Cat. No. LY-034

Description

Lysosomes are spherical vesicles in eukaryotic cells responsible for waste removal. The digestive enzymes in lysosomes play a vital role in digesting excess or worn-out organelles, food particles, and engulfed viruses and bacteria. Lysosomes are relatively large organelles ranging in size from 0.1 to 1.2 μm . The ability to isolate lysosomes is an essential first step in studying autophagy, protein degradation, and recycling in a cell. Traditional methods for isolating lysosomes are based on techniques developed in the 1970s, requiring a large amount of starting material and a tedious process. Significant cross-contaminations are inevitable. This kit employs a patented spin-column-based technology that is simple, rapid, and efficient. Native lysosomes can be significantly enriched without using a Dounce homogenizer and ultracentrifugation. The protocol takes less than 90 minutes, using only 20-30mg of the starting material.

Kit Components (20 prep)

1. Buffer A	15 ml
2. Buffer B	2 ml
3. Plastic rods	2
4. Filter Cartridge	20
5. Collection Tube	20

Additional Materials Required

1 X PBS, Vortexer, Table-Top Microcentrifuge with a maximum speed of 16,000 X g. **The centrifuge should be able to reach top speed within 10 seconds.**

Shipping and Storage: Ship at ambient temperature and store at 4°C

Important Information:

1. Read the entire protocol carefully. Chill protein extraction filter cartridge with collection tube on ice before use.
2. **All centrifugation steps should be performed at 4°C in a cold room or in a refrigerated microfuge.**
3. For protein phosphorylation studies, **phosphatase inhibitors** (such as PhosStop from Roche) should be added to buffer A prior to use. Add protease inhibitor cocktails to buffer A and B if protein degradation is a concern.



4. BCA Protein Assay Kit (Pierce, Cat #:23227) is recommended to determine protein concentration.

Protocol

1. Place the filter cartridges in the collection tubes and incubate them on ice.
2. **For cultured cells**, collect 25-30 X 10⁶ cells by low-speed centrifugation (500-600 X g for 5 min). Wash cells once with cold PBS, resuspend the pellet in 500 µl buffer A and incubate on ice for 5-10 min. **Vortex the tube vigorously for 10-30 seconds**, and immediately transfer to the filter cartridge. Go to step 4.
2. **For tissue samples**, place 20-30 mg of tissue (fresh or frozen) in a filter cartridge. Add 200 µl buffer A to the filter and grind the tissue with a plastic rod for 1 min by pushing the tissue against the surface of the filter repeatedly with twisting force. After grinding, add 300 µl buffer A to the same filter cartridge, pipette up and down to mix, and incubate on ice **with cap open** for 5 min. Go to step 4.
****The plastic rod is reusable. Clean with 70% alcohol or water.**
4. Cap the filter cartridge, invert a few times, and centrifuge at 16,000 X g for 30 seconds. (Optional: the flow through in the collection tube can be resuspended and re-passed through the same filter. This may increase the yield).
5. Discard the filter and resuspend the pellet by vigorous vortexing for 10 seconds. Centrifuge at 2000 X g for 3 min (the pellet contains nuclei, large cell debris, and un-ruptured cells if any).
6. Transfer all supernatant to a fresh 1.5 ml microfuge tube and centrifuge at 4°C for 15 min at 11,000 X g. The pellet contains mainly mitochondria and cell debris. After centrifugation, carefully transfer 400 µl supernatant to a fresh 1.5 ml tube and spin at 16,000 X g at 4°C for 30 min. Remove the supernatant completely.
7. Resuspend the pellet in 200 µl cold buffer A by pipetting up and down 60-100 times and vortex vigorously for 20 seconds. Centrifuge at 2000 X g for 4 min. Transfer the supernatant to a fresh 1.5 ml tube, add 100 µl buffer B, and vortex briefly to mix well (the supernatant to buffer B ratio is 2:1). Incubate the tube on ice for 30 min and centrifuge at 11,000 X g for 10 min. Remove all supernatant, and spin the tube at 11,000 X g briefly to bring down residual buffer; remove it altogether.
8. Resuspend the pellet in 50-150 µl PBS or buffer of your choice. This is the highly enriched lysosome fraction.

Tech Notes:

The Typical protein yield is 50-100 µg/sample. Increase the starting material if the target protein concentration is low. If there is no visible pellet in step 8, check the supernatant and pellet in step



6 for the presence of lysosomes after 11,000 X g spin. If the majority of lysosomes are found in the pellet, reduce the centrifugal force from 11,000 X g to 8,000-10,000 X g in step 6.

Depending upon downstream applications, the insoluble lysosome fraction can be dissolved in the following reagents.

Following protein solubilization reagents are recommended.

Product Name	Cat. No.	Applications
Minute™ Denaturing Protein Solubilization Reagent	WA-009	SDS-PAGE electrophoresis and Western blotting, trypsin digestion, purification of proteins with biotin labeling or histidine labeling, etc.
Minute™ Non-Denatured Protein Solubilization Reagent	WA-010	ELISA, immunoprecipitation/Co-IP, enzymatic activity determination and other applications.
Minute™ Protein Solubilization Reagent for MS	WA-011	Trypsin digestion and subsequent mass spectrometry analysis.