

# HYALURONIDASE ASSAY KIT, with Positive Control Enzyme

Assay kit for  
Screening of Hyaluronidase inhibitors, and  
Quantification of Hyaluronidase activity in cells and other biological  
fluids.

Quantity: 96 Tests

Cat.#: Ra003-01-HAK

DECEMBER 2021

**For Research Use Only.  
Not for use in diagnostic or therapeutic procedures.**

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0)728080244  
F: +31 (0)728080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

## PLEASE READ BEFORE USING ASSAY KIT

- **Avoid repeated freezing and thawing of the reagents and Enzyme preparations.**
- **Vials may contain small quantities of material, hence ensure that they are centrifuged prior to opening.**
- **Please read the whole manual as many solutions require preparation immediately before use i.e. preparation during the assay itself. The assay protocol indicates at which step certain solutions should be prepared.**

## Table of Contents

I. Background	3
II. Principles of the assay	4
III. Precautions and Limitations	5
IV. Reagents Provided	6
V. Materials required but not Supplied	6
VI. Reagent Preparation	7
VII HYALURONIDASE Activity in other biological samples	
A. Preparation of Cell culture	8
i- Adherent cells	8
ii- Non-Adherent Cells	8
B. Tissue Extract	8
C. Biological fluids	9
VIII. HYALURONIDASE Assay Protocol	9
A. Enzyme reaction	9
B. Detection	9
IX. HYALURONIDASE Inhibitor Screening Protocol	10
A. Concurrent inhibition	10
B. Pre-Reaction Inhibition	10
X. Data Interpretation	11
XI. References	11
XII. Troubleshooting	11
XIII. Related Products Available From Razie	12
XIV. Appendix, Reagent composition	12
XV. Ordering information	12

## 1. INTENDED USE

The hyaluronidase assay kit is an assay for quantitative detection of hyaluronidase in cell culture supernatants, human plasma, biological fluids and tissue samples. The hyaluronidase assay kit is for research use only. Not for use in diagnostic or therapeutic procedures.

### I. Background


Hyaluronidase are a group of enzymes that degrade hayluronic acid (HA). Also known as hyaluronan is composed of alternating residues of  $\beta$ -D-(1-3) glucuronic acid and  $\beta$ -D-(1-4)-N-Acetylglycosamine. HA is widely distributed in the body and forms one of the most versatile macromolecules of the extracellular matrix. Hyaluronan plays both structural and biological function through its ability to interact with other matrix components. The mammalian hyaluronidase cleaves hyaluronan and potentially similar glycosaminoglycans by hydrolysis of the 1-4 bond between N-acetyl glucosamine and D-glucuronic.

Traditionally measurement of the activity of these type of enzymes have been difficult to perform and standardize

The hyaluronidase assay by Razie is

- Non-radioactive,
- Fast and easy to perform,
- Sensitive and specific to hyaluronidase,
- Uses a universal 96-well plate format ideal for inhibitor studies

## PRINCIPLES OF THE TEST



The decrease in the OD of well D compared to A is directly proportional to the hyaluronidase activity.

### III. Precautions and Limitations

1. **For Research Use Only. Not for use in diagnostic or therapeutic procedures.**
2. The physical, chemical, and toxicological properties of the chemicals and reagents contained in the Hyaluronidase Assay Kit may not yet have been fully investigated. Therefore, Razie Ltd recommends the followings:
  - All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes wash immediately with water. In case of any hypersensitivity reaction see a medical doctor. Reagents are intended for research use only and are not for use in diagnostic or therapeutic procedures.
  - Do not mix or substitute reagents with those from other lots or other sources.
  - Do not use Positive control enzyme as a reference enzyme in any other tests or procedure
  - Do not use kit reagents beyond expiration date on label.
  - Do not expose kit reagents to strong light during storage or incubation.
  - Do not pipette by mouth.
  - Do not eat or smoke in areas where kit reagents or samples are handled.
  - Avoid contact of skin or mucous membranes with kit reagents or specimens.
  - Rubber or disposable latex gloves should be worn while handling kit reagents or specimens.
  - Reagents might contain strong preservatives and may be toxic if ingested.
  - Avoid contact of substrate solutions with oxidizing agents and metal.
  - Avoid splashing or generation of aerosols.
  - In order to avoid microbial contamination or cross-contamination of reagents or specimens which may invalidate the test use disposable pipette tips and/or pipettes.
  - Use clean, dedicated reagent trays for dispensing the substrate reagents.
  - Exposure to acids will inactivate the signal generating reagents.
  - Glass-distilled water or deionized water must be used for reagent preparation.
  - Substrate solutions must be at room temperature prior to use.
  - Decontaminate and dispose specimens and all potentially contaminated materials as if they could contain infectious agents. The preferred method of decontamination is autoclaving for a minimum of 1 hour at 121.5°C.
  - Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 1.0 % sodium hypochlorite. Allow 30 minutes for effective decontamination. Liquid waste containing acid must be neutralized prior to the addition of sodium hypochlorite.

#### IV. Reagents Provided

Catalog #	Component	Quantity	Storage Temperature
FBHAL01-6K	Reaction Buffer	6ml	4°C
FBHAL02-5K	Base Buffer	10ml	4°C
FBHAL-03-100	Additive 1	100 µl	4°C
FBHAL-04 BE-100	Positive Control	120 µl	-20°C
FBHAL-05-PB300	Buffer for Control Enzyme	200 µl	4°C
FB001-05-20YA	Strep-HRP	10 µl	-20°C
FB001-06-10K	HRP Substrate	10ml	4°C
RA1003--01- HAP	Hyaluronic acid -plate	1	-20°C

#### V. Materials Required But Not Supplied

- PBS, pH 7.4
- Tween 20
- Glass-distilled or deionized water (preferentially filtered through 0.45µm filter)
  
- 0.12 M HCl or H<sub>2</sub>SO<sub>4</sub> for stopping the peroxidase reaction.
- Phenylmethyl Sulfonyl Fluoride (PMSF) and other protease inhibitors (optional, to protect tissue/cell extract)
- 10 µl to 1,000 µl adjustable single channel micropipettes with disposable tips
  - 50 µl to 300 µl adjustable multichannel micropipette with disposable tips
  - Multichannel micropipette reservoir
- Centrifuge tubes, beakers, flasks (necessary for preparation of reagents or for dilution of enzyme/buffers).
- Glass Homogenizer (optional, for tissue extract assays)
  
- Device for delivery of wash solution such as multichannel wash bottle or automatic wash system (optional)

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
 184 Park Drive, Milton Park  
 Abingdon OX14 4SE.  
 T: +44 (0) 1235 828 200  
 F: +44 (0) 1235 820 482

North America   
 1035 Cambridge Street,  
 Cambridge, MA 02141.  
 T: +1 (617) 945-5033 or  
 T: +1 (800) 987-0985  
 F: +1 (617) 945-8218

Europe   
 Berenkoog 41,  
 1822 BH Alkmaar,  
 Netherlands  
 T: +31 (0)72 8080244  
 F: +31 (0)72 8080142

Switzerland   
 Via Lisano 3,  
 (CP.683)  
 CH-6900  
 T: +41 (0) 91 604 55 22  
 F: +41 (0) 91 605 17 85

Germany   
 Bockenheimer Landstr. 17/19  
 60325 Frankfurt/Main  
 T: +49 (0) 69 779099  
 F: +49 (0) 69 13376880

- 37°C incubator
- 96-well plate reader capable of reading at 450 nm (650 nm as optional reference wavelength).
- Inhibitors or cells/tissue.

## VI. Reagent Preparation

- **Please note this is how the reagents are prepared. Please refer to the actual protocol for when to prepare the reagents. Bring all reagents to the room temperature.**

### 1. Reaction buffer

#### Note:

The Reaction Buffer should be prepared immediately before use and used within one hour of preparation. A solid precipitate might form in additive 1 due to refrigeration. Bring the solution to the room temperature and vortex the tube vigorously to dissolve the precipitate. Briefly centrifuge the tubes prior to opening.

Mix 6ml of the Reaction Buffer (Cat# FBHAL 01-6K) with an equal volume (6ml) of filtered distilled water in a clean tube and vortex.

Add **80µl Additive 1** solution (Cat# FB 001-02-30) to the above diluted base buffer and vortex. Keep at room temperature.

### 2. Base Buffer (FBHAL 02-5K)

Base Buffer is ready for use, bring to room temperature just before using. Used to

- Dilute the test samples, any inhibitors to be tested
- Prepare cell / tissue extracts.

### 3. Positive control Enzyme (optional)

**(vial FBHAL-04-BE-100)** Defrost and pulse spin the vial in a microfuge **for 1 to 2 seconds only** to collate all liquid droplets. The positive control enzyme is ready for use. This enzyme should be used within one hour.

At the start of setting up the plate, add 50µl of the positive control **enzyme (FBHAL 04-BE-100)** to a well followed by 50µl **enzyme buffer (FBHAL-05-PB300)** to the same well. Mix the solution very gently by pipetting up and down. Repeat the same for the second designated positive control well. **There is no need for addition of any other buffer or additives** to these wells. Please note that the positive control will not be used in calculating results, it only shows that the assay has worked.

### 4. Strep-HRP

Immediately before use, dilute 8µl of the supplied Strep-HRP (Cat# FB001-05-20) in 12 ml **PBS**. Use the diluted Strep-HRP within 30 min. A total volume of 100 µl/well of diluted Strep-HRP is required in the assay.

## 5. HRP substrate

The supplied Peroxidase substrate solution (Cat# FB001-06-10K) is ready for use.

**Bring to room temperature prior to use.**

## 6. Controls

**i. Negative control (must be included in each test):** A negative control without hyaluronidase, test sample or cell extract should be included. This will provide the baseline reading (i.e. no enzyme present, hence strongest colour is produced in these wells.)

Add a further 50µl of the reaction buffer to four wells (i.e. 100µl of the reaction buffer is present in each well)

**ii. Positive control (optional):** See section VI, 3 above. **There is no need to add any other buffer or additives** to these two wells.

## VII. HYALURONIDASE Activity in other biological samples

### A. Preparation of Cell culture samples

#### i- Adherent cells:

Before analysis, it is important to be able to transfer your cell culture to a **Serum-Free Media (SFM)** as follows:

- 1- Grow your cell culture to confluence.
- 2-Discard the culture media containing serum.
- 3-Wash the cell culture 4X by adding 1X PBS to the culture vessel, gently shaking the culture vessel and discarding the PBS.
- 4-Add SFM to the culture vessel and culture overnight in the incubator. The culture time in SFM will depend on the nature of the cells and should be determined independently.
- 5- Collect the SFM from the culture vessel. Use aliquots of this SFM for the enzyme assay. Follow steps in enzyme assay as described (section VIII)

#### ii- Non-adherent cells:

1. Centrifuge  $2 \times 10^6$  to  $1 \times 10^7$  non-adherent cells grown in a suitable SFM to form a pellet at 1500g for 3 min (3000 RPM in an ordinary bench top centrifuge).
2. Centrifuge the supernatant (at 1500g for 2min) to pellet any cellular contaminant.

3. Treat this final supernatant as your test sample and carry out the hyaluronidase assay as described in section VIII.

### **B. Tissue Extract**

The following procedure must be carried out in a safety hood / cabinet. Appropriate eye protection and surgical mask should be used.

1. Wash the tissue thoroughly using cold PBS. The tissue must be clean and as far as possible free from contaminants such as blood or blood clots.
2. Cut a small piece of the tissue and place it inside a clean glass mortar. Add a volume of 1x Base Buffer (at least 1ml of buffer per every mg tissue)
3. Grind and homogenize the tissue.
4. Transfer the homogenate into a test tube or microfuge tube and centrifuge at around 2000g for 2min.
5. Use 50µl of the clear supernatant and carry out the enzyme assay as described (section VIII). It might be necessary to repeat centrifugation to obtain a supernatant free from tissue debris. **Avoid repeated cycles of freeze / thaw of this supernatant.**

### **C. Biological Fluids**

There is no particular preparation for biological fluids, as long as it is clear and free from any contaminants such as particulate material. This can be achieved as described for tissue extracts by a simple centrifugation step. Biological fluids should be kept ice cold and used immediately or frozen (at -80°C if possible) until use. **Avoid repeated cycles of freeze/thaw as this will lead to inactivation of the enzyme.**

For blood, it is advisable to use only plasma and not serum or whole blood to avoid contamination by activated platelets products.

## **VIII. HYALURONIDASE Assay Protocol**

**Bring the reaction buffer to the room temperature.**

### **A. Enzyme reaction**

1. Rehydrate the microwell plate by adding 100µl of PBS to each well of the plate. Incubate at room temperature for a minimum of 30 min or overnight at 4°C. Cover the microwell plate with the lid or a Plate Cover.

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0)72 8080244  
F: +31 (0)72 8080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

2. Remove the PBS from the microwells. Take care not to scratch the surface of the microwells. **Wash the wells once, by adding and removing 200 µl Distilled Water to the wells.** Ensure that all the liquid is removed by tapping plate onto paper towels. Use the microwell plate immediately or place upside down on a wet absorbent paper. Don't allow wells to dry.
3. Add 50µl of the reaction buffer to all wells. (positive control wells will not receive this, See VI.3 above)
4. Add 50 µl of the reaction buffer to four wells as the negative control.
5. Add 50µl of your test samples to duplicate wells (you can allocate more wells to your sample if you wish)
6. Incubate the microwell plate (**ideally at 37° C**) for **one hour** on a plate shaker.
7. Wash the wells four times with approximately 200 µl 1X PBS/ Tween20 (0.1% V/V) per well and then twice with approximately 200 µl Distilled Water with thorough aspiration of microwell contents between washes. Take care not to scratch the surface of the microwells.

You don't have to have a plate washer to use this kit. Washing stage for all sections of the manual can be carried out as described below:

a- Prepare the washing solution PBS/Tween20 as described in the manual. Pour enough of this solution into a flat bottom dish (around 30X15 cm) The solution should be around 5-6 cm deep.

b-Empty the plate at the end off the incubation period.

c- Submerge the plate into this wash solution. Make sure all wells are filled up with the wash solution. shake the plate inside the solution.

d- Rapidly shake empty the plate into a sink. Tap the plate dry on a pad of kitchen towels. This is one wash cycle. Repeat the above procedure for each wash cycle needed.

Tap the plate dry on a pad of kitchen towels! Use the plate immediately after washing or place upside down on a wet absorbent paper. Do not allow wells to dry.

(use multi channel pipette for partly covered plate, see appendix below)

## B. Detection

**Blank the plate reader according to the manufacturer's instructions against the air.**

1. Add 100 µl of diluted Strep-HRP at room temperature (prepared in section VI.4) per well. Incubate at **37° C** for 50 minutes on a plate shaker.
2. Wash the microwell plate five times with approximately 200 µl 1X PBS / well with thorough aspiration of microwell contents between washes. Take care not to scratch the surface of the microwells. After each wash, empty wells and tap microwell plate on paper towels to remove excess liquid. Use the plate immediately after washing or place upside down on a wet absorbent paper for not longer than 15 minutes. Do not allow wells to dry.
3. This section is carried out at room temperature. Add 100 µl of **room temperature** peroxidase substrate (see step VI.5) per well. Shake the plate by **Gently Tapping** the side of the plate to get a uniform colour. The colour of the liquid should start to turn blue. Please see the video clip on the AMSBIO site which shows an example of development stage. The colours generated in this clip are an example only and your test might show different colour intensity. In this clip, Wells 1A&B, negative control and 1C&D are positive controls.
4. Making sure that there no bubbles on the surface of the liquid, read absorbance on of each microwell on a spectrophotometer using 650 nm as the reference wavelength, at intervals of 1 to 2 minutes. This step should be carried at room temperature. Please note, the rate of colour generation is very dependent on the ambient temperature. The higher the temperature the faster colour develops. The colour could develop very fast within minutes.
5. The substrate reaction should be stopped as soon as an OD of 0.6 – 0.7 is reached for the blank wells. Stop the enzyme reaction by quickly pipetting 100 µl of 0.12M HCl into each well, including the blank wells. It is important that the Stop Solution is spread quickly and uniformly throughout the microwells to completely inactivate the enzyme. The colour of the liquid in wells should have turned into yellow at this stage.

**(Make sure there are not any bubbles at the surface of liquid)**

If you found colour development too fast, **use FB001-06-10K, HRP substrate, straight out of the fridge at 4 degrees or lower and dont let it warm up to room temperature.**

6. Determine the absorbance of all wells, at 450nm, within 5 min post stopping. The colour post stopping is relatively stable for one hour.

## IX. HYALURONIDASE Inhibitor Screening Assay Protocol

It is important to make sure that the inhibitor preparation is as pure and homogeneous as possible, so that the any effect seen, is representative of the inhibitor and not other compounds. The inhibitor preparation should be checked in

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0)72 8080244  
F: +31 (0)72 8080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880

the absence of any enzyme in this assay, to ensure that it has no effect on the mechanism of the assay. The wells should generate full colour with or without the inhibitor.

Similarly, for your targeted enzyme, it is recommended to run a pilot study to determine the enzyme's activity behavior, (e.g. check dilution response), using Razie's assay prior to the inhibition studies.

Potential inhibitors can be assayed using two different procedures

#### **A. Concurrent inhibition**

The inhibitor can be added to the enzyme at the time of preparing the reaction mixture.

- 1- Prepare dilutions of the inhibitor in the hyaluronidase buffer (Cat# FBHAL02-5K) you might be able to use another buffer provided that it does not have extreme pH or detectable levels of heavy ions.
- 2- Prepare a dilution of the enzyme from your source in the hyaluronidase buffer (prepared in step VI.2)
- 3- Follow procedure as VIII. A.1-5 above.
- 4- Add appropriate volumes of enzyme preparation and inhibitor together to one well. Make sure the final volume in the well is 100µl. (The base 50µl reaction buffer, step VIII. A.3, can be adjusted to achieve this volume).
- 5- Follow the procedure as for hyaluronidase assay (Steps VIII. A.6, 7).
- 6- Detection: Follow VIII B1-5

#### **B. Pre reaction inhibition**

In this procedure the inhibitor is added prior to the reaction.

- 1- Mix aliquots of the enzyme and the inhibitor as required (For example a ratio of 1:1 v/v i.e. 100µl of inhibitor to 100µl enzyme or 2:1, 3:1 and so on), mix.
- 2- Incubate at 4°C for the desired period of time. This time may vary between a few minutes to overnight. Please note that the longer the incubation period the more chance there is of natural enzyme degeneration.
- 3- Use 50µl of the above enzyme/inhibitor mixture as your test sample and follow the normal assay procedure as in VIII.A 1-7
- 4- Detection: Follow VIII B1-5

### **X. Data Interpretation.**

Calculate the activity of your sample (in ng HA removed per minute) using the following formulae

First calculate R,

$$R = (OD / \text{Max OD}) \times 500$$

MaxOD= maximum reading possible from your plate reader from the negative control wells (section VI, 6i and VIII),A4).

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
184 Park Drive, Milton Park  
Abingdon OX14 4SE.  
T: +44 (0) 1235 828 200  
F: +44 (0) 1235 820 482

North America   
1035 Cambridge Street,  
Cambridge, MA 02141.  
T: +1 (617) 945-5033 or  
T: +1 (800) 987-0985  
F: +1 (617) 945-8218

Europe   
Berenkoog 41,  
1822 BH Alkmaar,  
Netherlands  
T: +31 (0)728080244  
F: +31 (0)728080142

Switzerland   
Via Lisano 3,  
(CP.683)  
CH-6900  
T: +41 (0) 91 604 55 22  
F: +41 (0) 91 605 17 85

Germany   
Bockenheimer Landstr. 17/19  
60325 Frankfurt/Main  
T: +49 (0) 69 779099  
F: +49 (0) 69 13376880



OD= reading obtained from your test well

Then substitute R in the following formulae to obtain the activity

**Activity\* = (500- R) / 60 (ng HA removed per minute)**

Specific activity is then calculated by dividing the activity by protein concentration present in your sample to give: nanograms of HA removed per minute per mg protein of active enzyme.

\* This figure is based on the amount of HA loaded onto the plate and assumes that all have been immobilized by the active sites in the well.

**XI. Reference List**

1. M.R. Bansal, K.C. Kanwar and G.s. Gupta, Immunoenzymic Studies on testicular hyaluronidase from rhesus Monkeys (*Macaca mulatta*), *J.Reprod. Fert.* (1982) 64:267-273
2. D Watson, Hyaluronidase, *British Journal of Anaesthesia* (1993) 71: 422-425
3. Kazuhiro Mio et al. Detecting Hyaluronidase and hyaluronidase inhibitors, in *Proteoglycan Protocols, Methods in Molecular Biology book series, vol(171), Humana Press (2001), ISSN 1064-3745*
4. Krupkova, Olga, et al. "Expression and activity of hyaluronidases HYAL-1, HYAL-2 and HYAL-3 in the human intervertebral disc." *European Spine Journal* (2019) Nov 22

**XII. Troubleshooting**

PROBLEM	CAUSE	SOLUTION
No colour reaction is seen.	Strep-HRP has not been added.	Make sure step VIII, B.1 is carried out.
	Not enough incubation time	Prolong the incubation of strep-HRP.
No detectable signal is generated.	Cold HRP substrate has been used	Make sure HRP substrate is at room temperature
	Insufficient amount of the hyaluronidase enzyme is present	Increase the concentration of protein in your sample enzyme
		Increase incubation time to 1.5 hours.
	Increase the volume and/or concentration of cell extract added to each well.	
Colour generated is above the working	Long incubation of HRP substrate prior to stopping the reaction	Stop the reaction much earlier (Step VIII, B 4,5)

AMSBIO | [www.amsbio.com](http://www.amsbio.com) | [info@amsbio.com](mailto:info@amsbio.com)

UK & Rest of the World   
 184 Park Drive, Milton Park  
 Abingdon OX14 4SE.  
 T: +44 (0) 1235 828 200  
 F: +44 (0) 1235 820 482

North America   
 1035 Cambridge Street,  
 Cambridge, MA 02141.  
 T: +1 (617) 945-5033 or  
 T: +1 (800) 987-0985  
 F: +1 (617) 945-8218

Europe   
 Berenkoog 41,  
 1822 BH Alkmaar,  
 Netherlands  
 T: +31 (0)72 8080244  
 F: +31 (0)72 8080142

Switzerland   
 Via Lisano 3,  
 (CP.683)  
 CH-6900  
 T: +41 (0) 91 604 55 22  
 F: +41 (0) 91 605 17 85

Germany   
 Bockenheimer Landstr. 17/19  
 60325 Frankfurt/Main  
 T: +49 (0) 69 779099  
 F: +49 (0) 69 13376880

range of the plate reader.		Increase the dilution of Strep-HRP (Step VI, 4) (ie 6µl in 12ml)
----------------------------	--	--

### XIII. Related Products Available From Razie Biotech Ltd

Catalog #	Description	Quantity
Ra1003-01-HAP	HA coated plate	1 plate
FBHAL 02-5K	Buffer for dilution	1 Vial

### XIV. Appendix

#### Reagent composition:

- I. **1X PBS (pH 7.4) (not supplied):** 7.5 mM Na<sub>2</sub>HPO<sub>4</sub>, 2.5 mM NaH<sub>2</sub>PO<sub>4</sub>, 145 mM NaCl.
- II. **2X Base Buffer:** Proprietary buffer solution containing active ions.
- III. **Strep-HRP** Streptavidin Peroxidase tagged proprietary solution.
- IV. **HYALURONIDASE Buffer:** Proprietary buffer solution.
- V. **HRP Substrate:** 3,3',5,5'-Tetramethylbenzidine (TMB).

#### Using part of the plate:

You can cover part of the plate by using readily available masking tape, often used in decoration, see below picture.

The wells not used are covered by the masking tape for later use. The lid can still be used during the test and for storage. The washing stage of a partly masked plate should be done by using a multi channel pipette.

The masked plate can be stored at 4°C for extended period of time.

