

MA-0113

## ATP Assay

(100 wells, Colorimetric, OD 570 nm or Fluorometric, Ex/Em = 535/587 nm, Store at -20°C)

### Background Information:

ATP is the primary energy currency of living systems. Virtually all energy dependent processes utilize the chemical energy stored in the phosphate bond of ATP. ATP is formed exclusively in mitochondria and a variety of genetic diseases affect ATP formation in the mitochondria. There are many commercially available ATP assays which detect femtomoles or less of ATP using luminescence, but these kits require specialized instrumentation and utilize luciferase which can be difficult to maintain in active form. BioVision's ATP Colorimetric and Fluorometric Assay kit is designed to be a robust, simple method which utilizes the phosphorylation of glycerol to generate a product that is easily quantified by colorimetric (OD 570 nm) or fluorometric (Ex/Em = 535/587 nm) methods. The assay can detect down to 50 pmol (1  $\mu$ M) of ATP in various samples. The kit provides sufficient reagents for 100 assays.

### Assay Principle:

- 1 – ATP in the sample is used to phosphorylate glycerol
- 2 – Glycerol phosphate is oxidized with the formation of hydrogen peroxide
- 3 – The hydrogen peroxide is used by peroxidase to oxidize ADHP to form intense color and fluorescence

### Assay Components:

Assay Buffer	25 ml	WM	MA-0113-A
ADHP Solution	0.2 ml	Red	MA-0113-B
Glycerol Kinase	1 vial	Blue	MA-0113-C
Glycerol Phosphate Oxidase/HRP	1 vial	Green	MA-0113-D
ATP Standard	1 $\mu$ mol	Yellow	MA-0113-E

### User Supplied Materials:

- Homogenizer or Polytron
- Deproteinizing materials (See PI-0102, PI-0103)

### Storage and Handling:

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.

**Assay Buffer:** Warm to room temperature before use. Store at -20°C or 4°C.

**ADHP Solution:** Ready to use as supplied. DMSO freezes just below room temperature. Bring to room temperature before use. Store at -20°C.

**Glycerol Kinase:** Dissolve in 220  $\mu$ l Assay Buffer. Aliquot and store at -20°C.

**Glycerol Phosphate Oxidase/HRP:** Dissolve in 220  $\mu$ l Assay Buffer. Aliquot and store at -20°C.

**ATP Standard:** Dissolve in 100  $\mu$ l of DI water giving a 10 mM solution. Keep on ice while in use. Store at -20°C.

### ATP Assay Protocol:

**1. Sample Preparation:** Cells ( $10^6$ ) or tissues (10 mg) are homogenized or otherwise disrupted quickly in 100  $\mu$ l of a protein precipitant (TCA, PCA) on ice then centrifuged and the clear supernatant transferred to another tube. See details on protein removal in PI-0102 or PI-0103 datasheet. Transfer 2-10  $\mu$ l of sample to a 96-well plate. Adjust all well volumes to 50  $\mu$ l Assay Buffer. To correct for endogenous background, prepare each sample in pairs with one of the pair used as a non-ATP control.

#### Note:

- ATP is labile. Use only very fresh samples or snap freeze and store samples in liquid N<sub>2</sub>.
- Tissues samples contain enzymes that consume ATP rapidly. Work efficiently, keep samples cold, denature the enzymes as soon as possible after obtaining samples.

**2. Standard Curve Preparation:** For an absorbance-based assay, transfer 10  $\mu$ l of the ATP Standard to 240  $\mu$ l of diH<sub>2</sub>O and mix, giving a 400  $\mu$ M solution. Transfer 0–5–10–15–20–25  $\mu$ l of the standard into a series of wells and adjust all well volumes to 50  $\mu$ l with Assay Buffer giving 0–2–4–6–8–10 nmol ATP per well. For a fluorescence-based assay, dilute the ATP Standard 10X more (10  $\mu$ l to 90  $\mu$ l) with diH<sub>2</sub>O and transfer 0–5–10–15–20–25  $\mu$ l of the standard to a series of wells giving 0–200–400–600–800–1000 pmoles ATP per well. Adjust all well volumes to 50  $\mu$ l with Assay Buffer.

**3. Reaction Mix:** Prepare sufficient reagent for the number of samples and standards to be run:

For each well, prepare 50  $\mu$ l of Reaction Mix:


	Absorbance Based Assay	Bkgd Control	Fluorescence Based Assay	Bkgd Control
Assay Buffer	44 $\mu$ l	46 $\mu$ l	45.8 $\mu$ l	47.8 $\mu$ l
ADHP Solution	2 $\mu$ l	2 $\mu$ l	0.2 $\mu$ l	0.2 $\mu$ l
Glycerol Kinase	2 $\mu$ l	-----	2 $\mu$ l	-----
Glycerol Phosphate Oxidase/HRP	2 $\mu$ l	2 $\mu$ l	2 $\mu$ l	2 $\mu$ l


Mix well. Add 50  $\mu$ l of the Reaction Mix to each well containing ATP Standards or test samples. Add 50  $\mu$ l Background Control mix to all paired controls.

**4. Measurement:** Monitor the reaction progress with either absorbance (OD 570 nm) or fluorescence (Excitation/Emission = 535/587 nm). When the signal from the standards is no longer changing, the reaction is finished. It should take less than 30 minutes to reach endpoint.

### 5. Typical Results:

nmole Standard	Standard Raw Values		Background Corrected Values	
	OD	OD	OD	OD
0	0.0441	0.0436	0	0
2	0.2762	0.2808	0.2321	0.2372
4	0.5430	0.5413	0.4989	0.4977
6	0.8250	0.8043	0.7809	0.7607
8	1.1078	1.0725	1.0637	1.0289
10	1.3154	1.3335	1.2713	1.2899





**6. Calculation:** Subtract 0 Standard from all other standards. Plot the Standard Curve. Determine the slope of the standard curve. This defines the OD/nmol. Subtract each background control from its paired sample. Apply the slope of the standard curve to the background corrected sample values to get the sample ATP content in each well. To get the ATP content of the original raw samples:  
Background corrected ATP content in the well / volume added to well = nmoles ATP /  $\mu$ l of sample added to well  
nmoles ATP /  $\mu$ l of sample added to well X total volume of supernatant from sample disruption/centrifugation = nmole ATP in original sample.  
nmole ATP in original sample / mg of tissue (or # of cells) = nmole ATP / mg tissue (/# of cells)

***FOR RESEARCH USE ONLY! Not to be used on humans.***