

Human Albumin (ALB) Protein

Catalogue No.: abx670201

Human Albumin (ALB) Protein is a Native protein.

Target:	Albumin (ALB)
Origin:	Human
Tested Applications:	ELISA
Host:	Human
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	Unconjugated
Form:	Lyophilized
Storage:	Store at -20 °C after reconstitution. Avoid repeated freeze/thaw cycles.
UniProt Primary AC:	P02768 (UniProt , ExPASy)
Gene Symbol:	ALB
KEGG:	hsa:213
String:	9606.ENSP00000295897
Buffer:	Prior to lyophilization: Ammonium bicarbonate
Note:	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.