



DiscoveryPak™ Stem Cell Fate Regulator Set II

A convenient set of eight small molecules useful for stem cell differentiation, self-renewal and reprogramming.


Product	Biological Activity	Cat. No.	Size
BIX 01294	G9a HMTase inhibitor	1678-5	5 mg
CHIR 99021	GSK-3 β inhibitor	1677-5	5 mg
PD184352	MEK inhibitor	1585-5	5 mg
PD-0325901	MEK inhibitor	1643-2	2 mg
RG 108	DNMT inhibitor	1679-10	10 mg
SB-431542	ALK5 inhibitor	1674-1	1 mg
SU-5402	FGFR inhibitor	1645-05	500 μ g
VPA	HDAC1 inhibitor	1647-200	200 mg




BIX 01294




CHIR99021




PD184352




PD 0325901




RG 108



SB-431542



SU-5402



VPA

Cat. No.

K853-8

1 set

\$955.00