# Caspase-2/Ich1 Inhibitor, Z-VDVAD-FMK

**CATALOG #**: 1073-100

STORAGE CONDITIONS: Store at -20°C.

SHELF LIFE: 6 months under proper storage conditions

MOLECULAR WEIGHT: 695.7

**SEQUENCE**: Z-Val-Asp(OMe)-Val-Ala-Asp(OMe)-FMK

(FMK, fluoromethyl ketone)

**PURITY:** >98% by HPLC analysis.

PACKAGE CONTENTS: 100 ul (2 mM) Z-VDVAD-FMK

#### DESCRIPTION:

A synthetic peptide that irreversibly inhibits caspase-2 and related caspase activity. The inhibitor is designed as a methyl ester to facilitate cell permeability. (**CAUTION:** If the intended use is on purified or recombinant enzymes, esterase should be added to generate free carboxyl groups).

#### RECOMMENDED USAGE:

We recommend using 1000X dilutions for inhibiting caspase-2 activity in Jurkat cell culture (e.g., Add 1  $\mu$ l to 1 ml of culture medium). The optimal doses may vary for different cells and culture conditions.

#### REFERENCES:

- 1. Thornberry, N.A., et al. (1997) J. Biol. Chem. 272:17907-17911.
- 2. Talanian, R.V., et al. (1997) J. Biol. Chem. 272:9677-9682.

#### FOR RESEARCH USE ONLY! Not to be used in humans.

#### **RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- · Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

# Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

#### Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

## Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

# Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

#### Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

# Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

### **Growth Factors and Cytokines**

Monoclonal and Polyclonal Antibodies