BioVision Rev.01/10 For research use only

Cyclosporin A, Fusarium solani

CATALOG #: 1522-100 100 mg

1522-1g 1 gram

NH NH

STRCUTURE:

HIN HIN OH

MOLECULAR FORMULA: $C_{62}H_{111}N_{11}O_{12}$

MOLECULAR WEIGHT: 1202.6

CAS NUMBER: 59865-13-3

PHYSICAL APPEARANCE: White Power

PURITY: >99%

SOLUBILITY: Soluble in ethanol and DMSO; slightly soluble in

water.

STORAGE CONDITIONS: Store at +4 °C and keep dry.

DESCRIPTION: A fungal metabolite possessing potent

immunosuppressive properties. It inhibits the T-cell receptor signal transduction pathway via the formation of cyclosporin A-cyclophilin complex that inhibits calcineurin (protein phosphatase 2B. Cyclosporin A and its derivatives (including Alisporivir) inhibit multiple

coronaviruses .

CAUTION: TOXIC. Avoid contact by all modes of exposure.

FOR RESEARCH USE ONLY! Not to be used in humans!

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies