


EZSolution™ Y-27632

ALTERNATE NAME *trans*-4-[(1R)-1-Aminoethyl]-N-4-pyridinylcyclohexanecarboxamide dihydrochloride

CATALOG #: 1784-5 5 mg



MOLECULAR FORMULA: C₁₄H₂₁N₃O₂·2HCl

MOLECULAR WEIGHT: 320.26

CAS NUMBER: 129830-38-2

APPEARANCE: Liquid

FORMULATION: 10 mM solution in H₂O

PURITY: ≥98% by HPLC

STORAGE: Store at -20°C. Protect from air and light

DESCRIPTION: A 10 mM (5 mg in 1550 μl) solution of Y-27632, Dihydrochloride (Cat. No. 1596-5) in deionized water.

HANDLING: Do not take internally. Wear gloves and mask when handling the product! Avoid contact by all modes of exposure.

FOR RESEARCH USE ONLY! Not to be used in humans.

RELATED PRODUCTS:

- BIO (Cat. No. 1673-1)
- CHIR CHIR99021 (Cat. No. 1677-5, 25)
- GSK-3 Inhibitor, TWS119 (Cat. No. 1655-2)
- PD-0325901 (Cat. No. 1643-2)
- Pyrintegrin (Cat. No. 1629-1,5)
- EZSolution™ Pyrintegrin (Cat. No. 1737-1)
- SB-431542 (Cat. No. 1674-1)
- Thiazovivin (Cat. No. 1681-1,5)
- EZSolution™ Thiazovivin (Cat. No. 1736-1)
- Y-27632, 2HCl (Cat. No. 1596-1,5,50)