



HRP Conjugated Anti-HA Tag Mouse Monoclonal Antibody (4F6)

Cat #: ABT2045

Size: 50µl/200µl/200µl×5

Product Information

	Product Name: HRP Conjugated Anti-HA Tag Mouse Monoclonal Antibody (4F6)		
	Applications: WB		Isotype: Mouse IgG
	Reactivity: Null		
REF	Catalog Number: ABT2045	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: Null
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles		Note: Null

Background: Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA-molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:1000-1:3000).

Storage Buffer: Liquid in PBS, pH 7.4 and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.