

HRP Conjugated Anti-His Tag Mouse Monoclonal Antibody (5C3)

Cat #: ABT2055

Size: 50µl/200µl/200µl×5

Product Information

	Product Name: HRP Conjugated Anti-His Tag Mouse Monoclonal Antibody (5C3)		
	Applications: WB		Isotype: Mouse IgG
	Reactivity: Null		
REF	Catalog Number: ABT2055	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: Null
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles		Note: Null

Background: A polyhistidine-tag is an amino acid motif in proteins that consists of at least five histidine (His) residues, often at the N- or C-terminus of the protein. It is also known as hexa histidine-tag, 6xHis-tag, and by the trademarked name His-tag. Polyhistidine-tags are often used for affinity purification of polyhistidine-tagged recombinant proteins expressed in Escherichia coli and other prokaryotic expression systems.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:3000-1:10000).

Storage Buffer: Liquid in PBS, pH 7.4 and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.