

Malaria (Plasmodium Falciparum Antigen) Rapid Test Kit

Catalogue No.: abx090709

Malaria (Plasmodium Falciparum Antigen) Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Plasmodium Falciparum Antigen.

Target: Plasmodium Falciparum Antigen

Reactivity: Human

Tested Applications: LF

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Storage: Store at 2-30 °C.

Buffer: Not applicable.

Sample Type: Whole blood.

Note: This product is for research use only.
Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.

For Reference Only