

Adenovirus Antigen and Rotavirus Antigen Rapid Test Kit

Catalogue No.: abx090727

Adenovirus Antigen and Rotavirus Antigen Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the combined detection of Adenovirus Antigen and Rotavirus Antigen .

Target: Adenovirus Antigen and Rotavirus Antigen

Reactivity: Human

Tested Applications: LF

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Storage: Store at 2-30 °C.

Buffer: Not applicable.

Sample Type: Stool.

Note: This product is for research use only.
Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.

For Reference Only