

FACING THE  
FUTURE  
RESPONSIBLY

 **BIONOVAS**<sup>®</sup>

## Certificate of Analysis

**Product Name:** Agarose Low EEO  
**Cat #:** AA0481-0500  
**CAS #:** 9012-36-6  
**Lot #:** LOCT14  
**Storage:** Store at Room Temperature.

Appearance	Fine, white homogeneous powder
Moisture content	≤10%
Sulfate content	≤0.15%
Ash	≤2.0%
Gel Strength	1300g/cm <sup>2</sup>
EEO, mr	0.09-0.13
Gelling Temperature	36 ± 1.5 ° C
Melting Temperature	87-89 ° C
RNase/DNase	None Detected

 **BIONOVAS**

Email: [service@bionovas.com](mailto:service@bionovas.com)