

FACING THE  
FUTURE  
RESPONSIBLY

 BIONOVAS®

## CERTIFICATE OF ANALYSIS

Product	Coomassie Brilliant Blue R-250
Grade	Ultra Pure
Product Code	AC1270
Formula	$C_{45}H_{44}N_3O_7S_2Na$
MW	825.99
CAS#	6104-59-2
Lot No	KMAY13

Test Items	Specifications	Results
Appearance	Dark-Purple-Brown Powder	
Loss on drying	≤5.0%	
Dye content	≥80%	
Solubility (10mg/ml of water)	Dark blue-black solution	

Storage: RT. Protect from moisture.

 BIONOVAS

Email: [service@bionovas.com](mailto:service@bionovas.com)