

FACING THE
FUTURE
RESPONSIBLY

BIONOVAS®

CERTIFICATE OF ANALYSIS

Product	NP-40 (octylphenoxy polyethoxy ethanol)
Grade	Reagent
Product Code	AN0490-0500
Formula	N/A
MW	603
CAS#	9016-45-9
Lot No	

Test Items	Specifications
Appearance	Viscous liquid
Solubility (5 % W/V in Water)	Pass
Color Pt-Co	≤70
Specific Gravity	1.052 to 1.062
Cloud point (1%, aq.)	52 -56°C
Density	1.06 g/cm ³

Storage: 18~25°C

BIONOVAS

Email: service@bionovas.com