


Deoxymikanolide

Chemical Properties

CAS No. : 23753-57-3
 Formula : C₁₅H₁₆O₅
 Molecular Weight : 276.29
 Appearance : no data available
 Storage : Powder: -20°C for 3 years
 In solvent: -80°C for 2 years


Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6194	18.0969	36.1939
5 mM	0.7239	3.6194	7.2388
10 mM	0.3619	1.8097	3.6194
50 mM	0.0724	0.3619	0.7239

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481