


Dihydromikanolide

Chemical Properties

CAS No. : 23758-04-5
 Formula : C₁₅H₁₆O₆
 Molecular Weight : 292.29
 Appearance : no data available
 Storage : Powder: -20°C for 3 years
 In solvent: -80°C for 2 years



Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4213	17.1063	34.2126
5 mM	0.6843	3.4213	6.8425
10 mM	0.3421	1.7106	3.4213
50 mM	0.0684	0.3421	0.6843

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481