


## CERTIFICATE OF ANALYSIS

Product                   Agarose LM  
Product Code           AA0520  
CAS#                     9012-36-6  
Batch No                MJUL12



Test Items	Specifications	Results
Appearance	White to off-white powder	White powder
Gel strength (1% Gel)	>200 g/cm2	300g/cm2
EEO,-Mr	<0.1	0.07
Gelling Point (1.5% Gel)	28 °C ± 2.0 °C	29.5 °C
Melting Point (1.5% Gel)	<65 °C	64.5 °C
Sulfate	≤ 0.1%	0.09%
Water	≤ 10%	8.5%
DNase	Not Detected	Conform
RNase	Not Detected	Conform

Storage:    RT