

Dengue Virus Type 2, Envelope Protein

Catalogue No.: abx061506

Dengue Type 2 Envelope Protein (EP) is a recombinant protein produced in *Drosophila*.

Target:	Dengue Virus Type 2, Envelope
Expression:	Recombinant
Tested Applications:	ELISA, LF
Host:	<i>Drosophila</i>
Immunogen:	(a.a. 281-675)
Form:	Liquid
Purity:	> 85% (SDS-PAGE)
Purification:	Purified by affinity chromatography.
Storage:	Store at -20°C. Avoid repeated freeze/thaw cycles.
Molecular Weight:	50 kDa
Tag:	His tag
Buffer:	Phosphate Buffered Saline, pH 7.4. Preservative: 0.09% Sodium Azide.
Concentration:	2.8 mg/ml
Note:	This product is for research use only. This product is shipped with dry ice. Not for human consumption, cosmetic, therapeutic or diagnostic use.